

Graham Health System Laboratory
210 West Walnut or 180 S Main St
Canton, IL 61520
Phone; 309-647-5240 or 309-647-0201

PEDIATRIC URINE COLLECTION BAG
Patient collection handout

1. Separate the infant's legs. Make sure perianal and pubic areas are clean and dry.
2. Peel off backing paper on bag, exposing hypoallergenic adhesive.
3. For female use: Turn bag upside down with U shaped foam on bottom. Press adhesive firmly to skin surrounding urethral and vaginal openings, but not rectum. Fold bag forward.

For male use: Position penis naturally through bag opening with excess of bag toward the feet. Press adhesive firmly to entire perimeter without puckering.

4. Reapply diaper, if required, directly over bag.
5. After infant has urinated, carefully remove bag and place entire bag into urine collection cup.
6. Return specimen to lab within one hour. If unable to return specimen within one hour, refrigerate. Specimen is stable for 24 hours refrigerated.